


FLECOTEC

FLECOBOX (I-Box)



The table is optional and is available in various models for the I-Box

FLECOBOX I-Box

THIS INDIVIDUALLY PREFABRICATED FLEXIBLE ISOLATOR WITH SECONDARY CONTAINMENT AND INTEGRATED GLOVES MEETS ALL THE REQUIREMENTS OF CONTAINMENT LEVEL OEB 5 FOR HIGHLY TOXIC PRODUCTS.

The design of the I-Box does not require a base frame but is given its shape by the PE support tubes on the inside.

As a result and due to its very compact design, the I-Box is preferably used in laboratories where it offers personnel an additional high level of safety for a wide range of work steps.

STANDARD SIZES

| | 250 | 500 | 650 | 800 |
|----------------------|-----|-----|-----|-----|
| Width A (mm) | 250 | 500 | 650 | 800 |
| Length B (mm) | 450 | 500 | 650 | 650 |
| Height C (mm) | 450 | 500 | 650 | 650 |
| Base outlet possible | – | ✓ | ✓ | ✓ |
| Number of gloves (n) | 2 | 2-4 | 2-4 | 2-6 |

(special sizes on request)



Rommelag Flex Pharma GmbH

Breisgauring 13
79427 Eschbach
Germany

T + 49 7632 828 553-0
F + 49 7632 828 553-9

mail.rfp@rommelag.com
www.rommelag.com

